

Detail & Accessories

DEXCOWIN®

Mobile Hand-held X-ray system



Remote Control Connector

Anti Skid Grip

Exposure Button

Battery Slot & Cover

Time & Mode Control Button

LCD Window

Handle

BACK-SCATTER SHIELD



- For Radiation Safety**
- Radiation Below Detectable Limit (0.06 mR/hr) @1m From Focal Spot
 - Leakage Far Below FDA Limit (25 mR/hr) @1m From Focal Spot
 - Radiationh 2.5 mR/hr @5cm From Accessible Case Surface

REMOTE CONTROL



CARRYING CASE



BATTERY / CHARGER



- Material**
Charge Condition
 Current : 1.0 A
 Voltage : 16.8 VDC
 Temperature : 0 ~ 40 °C
- Use condition**
 Max. current : 26 A
 Voltage range : 12(min) ~ 16.8(MAX)V
 Temperature : -20 ~ 40 °C
 Full charging & discharging : 300 times
 Recharging cycle : 1.5 ~ 2.0 hrs

MOBILE STAND

